



# The Medical Desktop - Prescribing




Prescribing

## Prescribing



1. Click **Prescribe** button.
  2. **Search** for the **Medication** and select from the list.
  3. Click **OK**.
- The **Prescription Item Details** screen is displayed.
4. **Change Provider** if recording medication prescribed elsewhere. Also select **Given Else** tick box.



5. Enter **Dose**.
6. Select **Form**.
7. Select **Route**.
8. Select **Frequency**.
9. Enter **For** days and **Quantity**.
10. Enter **Repeats**.
11. Enter **Directions for Use**.
12. Click **Save** to add the dose to **Saved Defaults** for future prescriptions.
13. Click **OK**.

# The Medical Desktop - Prescribing



## Printing Prescriptions

The TOP – **Regular Medications**: Show patient's medications

The BOTTOM – **Medications History**: Shows the prescribing history of the medications.

The screenshot shows the 'Medications' window with two main sections: 'Regular Medications' and 'Medications History'. The 'Regular Medications' table has one entry for Buscopan Tablets. The 'Medications History' table shows three entries, including Glucose/Magnesium/Potassium/Sodium Chloride infusion, Diazepam, and Buscopan Tablets. On the right, a vertical toolbar contains buttons for 'Repeat', 'Repeat All', 'New', 'Modify', 'Delete', 'Print', and 'Close'. Red arrows and callouts indicate: 1. The 'Print' button; 2. The 'New' button; 3. The 'Repeat' button; 4. The 'Allergies...' button; 5. The 'Medications History' table.


Start Date	Expiry Date	Prescriber	Medication/Fomulation	Dose	Qty	Rp
20/05/2020		Ent.Fred	Buscopan Tablets , Tablets 10mg	1 Tablets Oral 2 Times Daily for 30 days 100 RP:2 Interva...	100	2

Date	Script No	Prescriber	Medication/Fomulation	Dose	Qty	Rp	Prt
29/04/2020	317134	Ent.Fred	Glucose; Magnesium; Potassium; Sodium; Chlor...	1 Infusion for 1 days 1	1	0	No
06/05/2020	974406	Ent.Fred	Diazepam (Tablets) 5mg	1 Tablets 500	500	0	No
20/05/2020	683934	Ent.Fred	Buscopan Tablets , Tablets 10mg	1 Tablets Oral 2 Times Daily for 30 days 100 RP:2 Interval...	100	2	No

To print a prescription...

1. Click the **Print** button.
2. When creating a NEW prescription from the **Medications** screen, use the bottom **New** button.  
**Note:** The top **NEW** button is for adding on patients existing medications.
3. Use the **Repeat** buttons to create a Repeat Prescription.
4. **Allergies** – Selecting **No Allergy** will display in green).
5. When you click **Print**, all items prescribed that day print on the one Rx form. The **Prt** column will display 'YES' to indicate the script has been printed.



The **Medications** button on the **Medical Desktop**, **F9** or the **Medications** button will display the Medications screen. To Prescribe a new medication, use the **New** button on the bottom half.



The **Prescribe** button on the **Medical Desktop**, will display your prescribing database first to enable selection of the medication.